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

——Using CLA as an Example

1. Basic Principles


ABE8e is an optimized adenine base editing system that enables the precise conversion of A•T base pairs to G•C base pairs without inducing double-strand breaks (DSBs) in DNA. ABE8e mainly consists of the following three functional modules: **(1) TadA8e**: An evolved adenine deaminase that deaminates adenine (A) on the target DNA strand, converting it into inosine (I), which is interpreted as guanine (G) during DNA replication. **(2) Cas9 nickase (D10A)**: A mutated version of Cas9 that introduces a single-strand nick on the non-edited DNA strand, working in conjunction with sgRNA to ensure site-specific targeting. **(3) sgRNA**: A single-guide RNA that directs Cas9 to the target DNA sequence, enabling the deaminase to act on the desired adenine base within a defined editing window.

Our vector system allows tandem integration of multiple **tRNA + target + gRNA** units (Fig. C), enhancing mutation efficiency. Plant endonucleases RNase P and RNase Z cleave the 5' and 3' ends of tRNA (Fig. A), releasing gRNAs that guide Cas9/gRNA to cleave the target DNA at the specified site.

In our lab, we modified the ABE8e vector by replacing the kanamycin resistance marker with hygromycin and incorporating the cotton endogenous U6 promoters (U6-7) for gRNA transcription. The protocol below focuses on constructing a vector with two target sites under the U6-7 promoter using rice Ubiquitin2 as the Cas9 promotor.

2. Target Search and Primer Design

CLA Sequence and Targets: Use the CRISPR-P website (<http://cbi.hzau.edu.cn/crispr/>) to search for targets (limited to *Gossypium raimondii* data). Alternatively, manually BLAST sequences for specificity. Prefer targets with unique restriction sites near the PAM (3 bp upstream) for downstream validation. T-to-C mutation on the negative strand corresponds to an A-to-G mutation on the positive strand. The target base is ideally located within positions 4 to 8 distal from the PAM site.

CTCGACGGACCTATACCACCTGTGGGAGCTTTGAGCAGTGCTCTCAGCAGGCTGCAAT
CAAACAGGCCTCTTAGAGAACTGAGAGAGGTTGCAAAGGTAAGTAGAAATGAATAGA
ACCCAGAAACTAATATAATTTTACTAGAACTAAACTTACTAAGACTGGATGTGGT
TTGATTGCAGGGAGTTA **CCAAGCAATCGGTGGGCC** **TATGCATGA**ACTGGCTG

CAAAAGTTGATGAGTATGCTCGTGAATGATCAGTGTTCTGGATCTACACTTTTCGAA

GAACCTGGACTATATTATATTGGACCTGTTGATGGCCACAACATCGATGATTTAGTTTCT
 ATTCTCAAAGAGGTTAAGACTACTAAAACAACGGGTCCGGTCTTGATCCATGTTGTCAC
 TGAGAAAGGCCGAGGTTATCCATATGCAGAGAGAGCTGCTGACAAGTATCACGGTAAC
 ATACAGAATAAGCCTTTTTAGAAGAGAGATCATTTTCATATGATTTAATTTACTGGTGCCT
 CGATATCTGACCTTAGTTTAGGGAATAAAAGAATAAACTGTACCTGGATTCTTTTCTCA
 GGAGTGGTG**AGTTTCGATCCGGCAAC**TGGAAAGCAATTCAAAGGCAA

TTCTGCTACCCAGTCTTACACTACATATTTTGCTGAGGCTTTGATTGCGGAAGCTGAGG
 CAGACAAAAATATTGTTGCCATCCATGCTGCAATGGGAGGTGGAACCGGATTAAACCT
 CTTCTCCGCCGGTTCCCAAAAGATGTTTCGATGTGGGGATAGCTGAACAACATGCTG
 TCACCTTTGCTGCAGGCTTGGCCTGTGAAGGCTTGAAACCTTTTTGTGCAATCTACTCA
 TCATTCATGCAAAGGGCTTATGACCAG

Insert the target sequence and ligate between the MfeI digestion sites shown in the box of the vector diagram. The final ligation method is the one-step cloning method.

*aagcatcagatgggca***AACAAAGCACCAGTGGTCTAGTGGT**AGATAAGTACCCTGC
CACGGTACAGACCCGGGTTTCGATTCCCGGCTGGTGCA**CATAGGCCACCGA**
TTTGCTGTTTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTA
TCAACTTGAAAAAGTGGCACCAGTCGGTGC**AACAAAGCACCAGTGGT**
CTAGTGGTAGATAAGTACCCTGCCACGGTACAGACCCGGGTTTCGATTC**CCG**
GCTGGTGCA**GTGAAGTTCGATCCGGCAAC***gtttagagctagaa*..... (载体 gRNA)

notes: Note: Lowercase = vector backbone; underscores = tRNA; bold = gRNA; italics = target site.

Primers::

CLA1 as: 5' AGCAAATCGGTGGGCCTATGtcaccagccgggaat3' (靶标 1 互补序列)
 CLA 2 s: 5' CATAGGCCACCGATTTGCTgttttagagctagaata 3' (靶标 1 序列)
 CLA 2 as: 5' GTTGCCGGATCGAACTTCACtgcaccagccgggaat3' (靶标 2 互补序列)
 inf CLA as: 5' ttctagctctaaaacGTTGCCGGATCGAACTTCAC3' (靶标 2 互补序列)
 pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAAGCACCAAGTGGTCTAG3'
 inf pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAA 3'



The lowercase part is a fixed sequence. Primer design only requires modifying the target sequence indicated by the underline in the primer.

Among them, inf pRGE32-7 s, pRGE32-7 s, u6-7 s are universal primers and do not need to be changed. pRGE32-7s contain a part of the vector sequence, facilitating In-fusion ligation. Because there are multiple trnas in the inserted sequence, inf pRGE32-7 s, in order to enhance amplification specificity, have shorter sequence lengths than pRGE32-7 s. Inf CLA as contains carrier sequence connectors for convenient subsequent In-fusion connections.

3. First PCR Amplification

Overlap Extension PCR:

This method generates fusion fragments by designing primers with complementary overhangs. Two overlapping PCR products are combined and extended to form a single construct.


The gRNA+tRNA combination sequence can be synthesized by itself and ligated to a common vector (such as using the T vector). In the first PCR: Fragments 1 and 2 were amplified respectively from the bacterial liquid or plasmid containing the gRNA+tRNA vector. (The relevant sequence can be found in the attachment.)

Reaction Mixture (20 μ L):

	1	2
ddH ₂ O	16.1	16.1
Buffer	2	2
dNTP	0.3	0.3
S primer	pRGEb32-7 s 0.2	CLA 2 s 0.2
AS primer	CLA1 as 0.2	CLA 2 as 0.2
Taq	0.2	0.2
Template DNA	1	1

PCR Conditions:

1. Initial denaturation: 95° C for 4 min
2. Denaturation: 95° C for 30 s
3. Annealing: 55° C for 30 s
4. Extension: 72° C for 20 s
 - Repeat steps 2 – 4 for 3 cycles
5. Denaturation: 95° C for 30 s
6. Annealing: 60° C for 30 s
7. Extension: 72° C for 20 s
 - Repeat steps 5 – 7 for 27 cycles
8. Final extension: 72° C for 5 min

Electrophoresis: Verify amplification.


4. Second PCR and Purification

Combine the two PCR products from Step 3 using overlap extension,

Reaction Mixture (100 µL)			
ddH ₂ O			83.5µl
Buffer			10µl
dNTP			1.5µl
S primer	Inf	pRGEB32-7 s	1µl
AS primer	Inf	CLA as	1µl
Taq			1µl
Fragment 1			1µl
Fragment 2			1µl

PCR Conditions:

1. Initial denaturation: 95° C for 4 min
2. Denaturation: 95° C for 30 s

3. Annealing: 59° C for 30 s
4. Extension: 72° C for 20 s
 - Repeat steps 2 – 4 for 28 cycles
5. Final extension: 72° C for 5 min

Purification: Use a gel extraction kit to purify the product and measure concentration.

5. Vector Digestion

Digest the vector with BsaI.

Reaction Mixture (100 μ L):

	Volume (μ L):
ABE8e	10 μ g
CutSmart Buffer	10 μ L
MfeI	1 μ L
ddH ₂ O	Up to 100 μ L

Digestion: 37°C for 5.5 h.

6. In-Fusion Cloning

Assemble the digested vector and insert using In-Fusion technology.

Components:

- **Insert (100 ng)**
- **Linearized vector (100 ng)**
- **Exnase (0.5 μ L)**
- **CE Buffer (1 μ L)**

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

7. Transformation

Transform competent cells with the ligation product.

Steps:

Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7s: TGTGCCACTCCAAAGACATCAG + Inf XX AS


Positive clones are sequenced (GC-rich regions).

Electrophoresis: Confirm positive clones.

Positive detection electrophoresis pattern:

Construct a successful carrier diagram:

8. Agrobacterium Electroporation

Electroporate the verified construct into *Agrobacterium*.

u6-7 s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

9. Validation

Genomic DNA Extraction: PCR amplify the target region, clone, and sequence.