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


——Using CLA as an Example

## 1. Basic Principles


After the target sequence is cleaved by the CRISPR system, the organism initiates its own repair mechanisms, leading to mutations and silencing of the target gene function (as shown below). The Cas9/gRNA complex can target any DNA site containing the motifs 5'-N<sub>20</sub>-NGG-3' or 5'-CCN-N<sub>20</sub>-3' (where NGG is the PAM required for Cas9 recognition). Cas9/gRNA precisely cleaves the target DNA double strand approximately 3 bp upstream of the PAM.

Our vector system allows tandem integration of multiple **tRNA + target + gRNA** units (Fig. C), enhancing mutation efficiency. Plant endonucleases RNase P and RNase Z cleave the 5' and 3' ends of tRNA (Fig. A), releasing gRNAs that guide Cas9/gRNA to cleave the target DNA at the specified site.

In our lab, we modified the pRGEB32 vector by replacing the kanamycin resistance marker with hygromycin and incorporating the cotton endogenous U6 promoters (U6-7, U6-9) for gRNA transcription. The protocol below focuses on constructing a vector with two target sites under the U6-7 promoter using rice Ubiquitin2 as the Cas9 promoter.

## 2. Target Search and Primer Design

**CLA Sequence and Targets:** Use the CRISPR-P website (<http://cbi.hzau.edu.cn/crispr/>) to search for targets (limited to *Gossypium raimondii* data). Alternatively, manually BLAST sequences for specificity. Prefer targets with unique restriction sites near the PAM (3 bp upstream) for downstream validation.

CTCGACGGACCTATACCACCTGTGGGAGCTTTGAGCAGTGCTCTCAGCAGGCTGCAAT  
CAAACAGGCCTCTTAGAGAACTGAGAGAGGTTGCAAAGGTAAGTAGAAATGAATAGA  
ACCCAGAAACTAATATAATTTTACTAGAACTAAACTTACTAAGACTGGATGTGGT  
TTGATTGCAGGGAGTTA **CCA**AGCAAATCGGTGGGCCTATGCATGAACTGGCTG

CAAAAGTTGATGAGTATGCTCGTGGAATGATCAGTGGTTCTGGATCTACACTTTTCGAA  
 GAACTTGGACTATATTATATTGGACCTGTTGATGGCCACAACATCGATGATTTAGTTTCT  
 ATTCTCAAAGAGGTTAAGACTACTAAAACAACGGGTCCGGTCTTGATCCATGTTGTCAC  
 TGAGAAAGGCCGAGGTTATCCATATGCAGAGAGAGCTGCTGACAAGTATCACGGTAAC  
 ATACAGAATAAGCCTTTTTAGAAAGAGAGATCATTTTCATATGATTTAATTTACTGGTGCCT  
 CGATATCTGACCTTAGTTTAGGGAATAAAAGAATAAACTGTACCTGGATTTCTTTTCTCA  
 GGAGTGGTGAAGTTCGATCCGGCAAC **TGG**AAAGCAATTCAAAGGCAA

TTCTGCTACCCAGTCTTACACTACATATTTTGCTGAGGCTTTGATTGCGGAAGCTGAGG  
 CAGACAAAAATATTGTTGCCATCCATGCTGCAATGGGAGGTGGAACCGGATTAAACCT  
 CTTCTCCGCCGGTCCCACAAAGATGTTTCGATGTGGGGATAGCTGAACAACATGCTG  
 TCACCTTTGCTGCAGGCTTGGCCTGTGAAGGCTTGAAACCTTTTGTGCAATCTACTCA  
 TCATTCATGCAAAGGGCTTATGACCAG

Insert the target sequence and ligate between the two Bsa I digestion sites shown in the box of the vector diagram. The final ligation method is the one-step cloning method.

*aagcatcagatgggca* **AACAAAGCACCAGTGGTCTAGTGGTAGAATAGTACCCTGC**  
**CACGGTACAGACCCGGGTTCGATTCCCGGCTGGTGCA** *CATAGGCCACCGA*  
*TTTGCT* **GTTTTAGAGCTAGAAATAGCAAGTTAAATAAGGCTAGTCCGTTA**  
**TCAACTTGAAAAAGTGGCACCAGTCCGGTGC** **AACAAAGCACCAGTGGT**  
**CTAGTGGTAGAATAGTACCCTGCCACGGTACAGACCCGGGTTCGATTCCCG**  
**GCTGGTGCA** *GTGAAGTTCGATCCGGCAAC* *gttttagagctagaa*..... (载体 gRNA)

notes: Note: Lowercase = vector backbone; underscores = tRNA; bold = gRNA; italics = target site.

#### Primers::

CLA1 as: 5' AGCAAATCGGTGGGCCTATG *tcaccagccgggaat3'* (靶标 1 互补序列)  
 CLA 2 s: 5' CATAGGCCACCGATTTGCT *gttttagagctagaaata 3'* (靶标 1 序列)  
 CLA 2 as: 5' GTTGCCGGATCGAACTTCAC *tcaccagccgggaat3'* (靶标 2 互补序列)  
 inf CLA as: 5' *ttctagctctaaac* GTTGCCGGATCGAACTTCAC *3'* (靶标 2 互补序列)  
 pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAAGCACCAGTGGTCTAG3'  
 inf pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAA 3'

#### (U6-9

pRGE32-9 s: 5' CAGCACATAACTGGCAAACAAAGCACCAGTGGTCTAG3'  
 inf pRGE32-9 s: 5' CAGCACATAACTGGCAAACAAA 3')

**The lowercase part is a fixed sequence. Primer design only requires modifying the target sequence indicated by the underline in the primer.**


Among them, inf pRGE32-7 s, pRGE32-7 s, pRGE32-9 s, inf pRGE32-9 s, u6-7 s, and u6-9 s are universal primers and do not need to be changed. pRGE32-7s and pRGE32-9s contain a part of the vector sequence, facilitating In-fusion ligation. Because there are multiple trnas in the inserted sequence, inf pRGE32-7 s and inf pRGE32-9 s, in order to enhance amplification specificity, have shorter sequence

lengths than pRGEB32-7 s and pRGEB32-9 s. Inf CLA as contains carrier sequence connectors for convenient subsequent In-fusion connections.


### 3. First PCR Amplification

#### Overlap Extension PCR:

This method generates fusion fragments by designing primers with complementary overhangs. Two overlapping PCR products are combined and extended to form a single construct.

Overlapping extension schematic diagram

This is a schematic diagram of the test process

The gRNA+tRNA combination sequence can be synthesized by itself and ligated to a


common vector (such as using the T vector). In the first PCR: Fragments 1 and 2 were amplified respectively from the bacterial liquid or plasmid containing the gRNA+tRNA vector. (The relevant sequence can be found in the attachment.)

**Reaction Mixture (20  $\mu$ L):**

|                         | 1                      | 2                   |
|-------------------------|------------------------|---------------------|
| <b>ddH<sub>2</sub>O</b> | <b>16.1</b>            | <b>16.1</b>         |
| <b>Buffer</b>           | <b>2</b>               | <b>2</b>            |
| <b>dNTP</b>             | <b>0.3</b>             | <b>0.3</b>          |
| <b>S primer</b>         | <b>pRGEB32-7 s 0.2</b> | <b>CLA 2 s 0.2</b>  |
| <b>AS primer</b>        | <b>CLA1 as 0.2</b>     | <b>CLA 2 as 0.2</b> |
| <b>Taq</b>              | <b>0.2</b>             | <b>0.2</b>          |
| <b>Template DNA</b>     | <b>1</b>               | <b>1</b>            |

**PCR Conditions:**

1. Initial denaturation: 95°C for 4 min
2. Denaturation: 95°C for 30 s
3. Annealing: 55°C for 30 s
4. Extension: 72°C for 20 s
  - Repeat steps 2–4 for 3 cycles
5. Denaturation: 95°C for 30 s
6. Annealing: 60°C for 30 s
7. Extension: 72°C for 20 s
  - Repeat steps 5–7 for 27 cycles
8. Final extension: 72°C for 5 min

**Electrophoresis:** Verify amplification.


#### 4. Second PCR and Purification

Combine the two PCR products from Step 3 using overlap extension,

| Reaction Mixture (100 $\mu$ L) |     |            |              |
|--------------------------------|-----|------------|--------------|
| ddH <sub>2</sub> O             |     |            | 83.5 $\mu$ l |
| Buffer                         |     |            | 10 $\mu$ l   |
| dNTP                           |     |            | 1.5 $\mu$ l  |
| S primer                       | Inf | pRGE32-7 s | 1 $\mu$ l    |
| AS primer                      | Inf | CLA as     | 1 $\mu$ l    |
| Taq                            |     |            | 1 $\mu$ l    |
| Fragment 1                     |     |            | 1 $\mu$ l    |
| Fragment 2                     |     |            | 1 $\mu$ l    |

#### PCR Conditions:

1. Initial denaturation: 95°C for 4 min
2. Denaturation: 95°C for 30 s
3. Annealing: 59°C for 30 s
4. Extension: 72°C for 20 s
  - Repeat steps 2–4 for 28 cycles
5. Final extension: 72°C for 5 min

**Purification:** Use a gel extraction kit to purify the product and measure concentration.

#### 5. Vector Digestion

**Digest the vector with BsaI.**

**Reaction Mixture (100  $\mu$ L):**

| <b>Volume (<math>\mu</math>L):</b> |                   |
|------------------------------------|-------------------|
| PCR Product                        | 10 $\mu$ g        |
| CutSmart Buffer                    | 10 $\mu$ L        |
| BSA                                | 1 $\mu$ L         |
| ddH <sub>2</sub> O                 | Up to 100 $\mu$ L |

Digestion: 37°C for 5.5 h.

**6. In-Fusion Cloning**

Assemble the digested vector and insert using In-Fusion technology.

**Components:**

- Insert (100 ng)
- Linearized vector (100 ng)
- Exnase (0.5  $\mu$ L)
- CE Buffer (1  $\mu$ L)

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

**7. Transformation**

Transform competent cells with the ligation product.

**Steps:**

Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

u6-9s: GTCAAAAACCTATCCCACATTGCTAA + Inf XX AS


Positive clones are sequenced (GC-rich regions).

Electrophoresis: Confirm positive clones.

Positive detection electrophoresis pattern:

### Construct a successful carrier diagram:

### 8. Agrobacterium Electroporation


Electroporate the verified construct into *Agrobacterium*.

u6-7 s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

u6-9 s: GTCAAAACTATCCCACATTGCTAA + Inf XX AS

### 9. Validation

1. Genomic DNA Extraction: PCR amplify the target region, clone, and sequence.
2. Restriction Enzyme Digestion: Validate via PAM-proximal restriction sites.

### 3.T7EI Assay: Detect mismatches in edited sequences.