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Design custom gRNAs Lookup human transcript expression Human Model organisms RNA viruses Additional resources

Design custom gRNAs

To accurately predict local target site accessibility use sequences of > 80 nt length.
To limit run time when designing guides on webpage, please use target sequences of < 500 nt length. For greater length requirements, please use the 'upload multiple single-entry fasta files' option.
Please be patient while the gRNAs are being scored. (100nt ~ 10sec; 500nt ~ 30sec)

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
pRGEB32-7 s: 5’ aagcatcagatgggcaAACAAAGCACCAGTGGTCTAG3’

Inf-cla as: 5’ atttgtgacctgcaCTGTAATCACCCCACAAATC 3'

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100μl 体系:	
Cas13	10μg
10X cut smart Buffer	10μl
BSA1	2 μl
Sbf1	2 μl
ddH ₂ O	up to 100μl

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u6-7 s: TGTGCCACTCCAAAGACATCAG + Inf XX AS


u6-9 s: GTCAAAAACACTATCCCACATTGCTAA + Inf XX AS

CRISPR/Cas13 protocol (英文版)


——Using CLA as an Example

1. Basic Principles


Cas13/crRNA can target any transcript/RNA for degradation.

Our vector system allows tandem integration of multiple **tRNA + target + crRNA** units (Fig. C), enhancing mutation efficiency. Plant endonucleases RNase P and RNase Z cleave the 5' and 3' ends of tRNA (Fig. A), releasing crRNAs that guide Cas13/crRNA to cleave the target RNA at the specified site.

In our lab, we modified the pRGEB32 vector by replacing the kanamycin resistance marker with hygromycin and incorporating the cotton endogenous U6 promoters (U6-7) for crRNA transcription. The protocol below focuses on constructing a vector with two target sites under the U6-7 promoter using rice Ubiquitin2 as the Cas13 promotor.

2. Target Search and Primer Design

CLA Sequence and Targets: Use the Cas13 design website (<https://cas13design.nygenome.org/>) to search for targets . Alternatively, manually BLAST sequences for specificity.

cas13design: A flexible tool to design Cas13d guide RNAs

 Design custom gRNAs

 Lookup human transcript expression

 Human

 Model organisms

 RNA viruses

 Additional resources

Design custom gRNAs

To accurately predict local target site accessibility use sequences of > 80 nt length.

To limit run time when designing guides on webpage, please use target sequences of < 500 nt length. For greater length requirements, please use the 'upload multiple single-entry fasta files' option.

Please be patient while the gRNAs are being scored. (100nt ~10sec; 500nt ~ 30sec)

Select Input Type:

- ☒ Paste single target
- ☐ Upload one single-entry fasta file
- ☐ Upload multiple single-entry fasta files (results e-mailed)

Enter target RNA sequence 

```
GCGGCTCGGGACGGAGGACGCGCTAGTGTCTTCTGT  
GTGGCAGTT CAGAAATGATGGATCAAGCTAGATCAGCATT  
CTCTAACTTGTGGTGGAACC
```

 Design guides

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
Cas13 guide RNA predictions

Please use the search field for searching or rank the guide RNA predictions column-wise.

 Download predictions

ATGGCCCTTTGTGCATCTTCATTTCTGCGCATTATCAACTGGGGTGCAGCTT
 CAGATCCTCAAAAATCCACTCCTTTTGCTTCCCATTTTCTTGGTGGATCAGA
 TCTGGTTTTGCAGCCATTGAAGAAGCTCAATCAGGTCAAGAAAAGGCCAG
 GTGGGGTTTATGCATCACTATCAGAAGGAGCTGAATATCACTCCCAAAGAC
 CAGCAACGCCTCTCTTGGACACTATAAACTATCCAATTCATATGAAAAATCT
 CTCTGTCAAGGAACTGAAACAACCTATCTGAAGAACTGCGGTCTGATGTGAT
 TTTCAATGTTTTCAAAAACCTGGGGGTCACTTGGGTTCAAGCCTTGGTGTGGT
 TGAAGTCACTGTGGCTCTTCATTATGTCTTCAATGCCCTAGAGACAAGATA
 TTGTGGGATGTTGGTCATCAGTCTTACCCTCACAAAATCTTGACTGGGAGA
 AGACATAAGATGCATACCATGAGGCCAACTAATGGATTGGCCGGATTACG
 AAACGGTCGGAGAGTGAATACGATTGCTTCGGGACTGGTCACAGTTCAAC
 CACAATCTCTGCAGGCTTGGGGATGGCTGTGGGAAGGGATCTGAAAGGTG
 AAAGGAATCATGTTGTTGCTGTCATAGGTGATGGTGCTATGACTGCAGGAC
 AAGCTTACGAAGCAATGAACAATGCCGGATATCTTGATTCTGATATGATTGT
 TATTCTTAATGACAATAAACAAGTTTCTCTGCCAACTGCCAATCTC

Insert the target sequence and ligate between the two Bsa I digestion sites and SbfI shown in the box of the vector diagram. The final ligation method is the one-step cloning method.

crRNA direct repeat/Cas13x.1

GCTGGAGCAGCCCCCGATTTGTGGGGTGATTACAGC

*aagcatcagatgggca***AACAAAGCACCAGTGGTCTAGTGGTAGAATAGTACCCTGC**
CACGGTACAGACCCGGGTTTCGATTCCCGGCTGGTGCA**CATAGGCCACCGA**
TTTGCTGCTGGAGCAGCCCCCGATTTGTGGGGTGATTACAGC**AACAAAG**
CACCAGTGGTCTAGTGGTAGAATAGTACCCTGCCACGGTACAGACCCGGGT
TCGATTCCCGGCTGGTGCA**GTGAAGTTCGATCCGGCAAC****GCTGGAGCAGCC**
CCCGATTTGTGGGGTGATTACAGC*tgagggtccacaaat.....*

notes: Note: Lowercase = vector backbone; underscores = tRNA; bold = crRNA; italics = target site; crRNA in red.

Direct synthesis is recommended, single-target reference sequences are available upon request.

*aagcatcagatgggca***AACAAAGCACCAAGTGGTCTAGTGGTAGAATAGTACCCTGC**
CACGGTACAGACCCGGGTTCGATTCCCGGCTGGTGCA*CATAGGCCACCGA*
*TTTGCT***GCTGGAGCAGCCCCGATTTGTGGGGTGATTACAGC***tgcaggccac*
aat.....

pRGE32-7 s: 5' *aagcatcagatgggca***AACAAAGCACCAAGTGGTCTAG**3'

Inf-cla as: 5' *atttgtagacctgca***CTGTAATCACCCACAAATC** 3'

3. Vector Digestion

Digest the vector with BsaI.

Reaction Mixture (100 μ L):

<i>Volume (μL):</i>	
Cas13	10 μ g
10X cut smart Buffer	10 μ l
BSA1	2 μ l
Sbf1	2 μ l
ddH ₂ O	up to 100 μ l

Digestion: 37°C for 5.5 h.

4. In-Fusion Cloning

Assemble the digested vector and insert using In-Fusion technology.

Components:

- Insert (100 ng)
- Linearized vector (100 ng)
- Exnase (0.5 μ L)
- CE Buffer (1 μ L)

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

5. Transformation

Transform competent cells with the ligation product.


Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7 s: 5' TGTGCCACTCCAAAGACATCAG 3'

Ubi-as: 5' TGTTGGTCGCCGTTAGGA 3'

Construct a successful carrier diagram:

6. Agrobacterium Electroporation

Electroporate the verified construct into *Agrobacterium*.