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
——Using CLA as an Example


1. Basic Principles

Tad-CBE is an optimized cytosine base editing system capable of precisely converting C•G base pairs into T•A base pairs without introducing DNA double-strand breaks (DSBs). Tad-CBE consists of the following three functional modules:


- (1) **Tad deaminase**: deaminates the target cytosine (C) on the DNA strand, converting it into thymine (T);
- (2) **Cas9 nickase (D10A)**: a mutated Cas9 enzyme that generates a single-strand nick on the non-edited strand, guided by the sgRNA;
- (3) **sgRNA**: guides Cas9 to accurately recognize the target sequence, enabling the deaminase to act on the cytosine base within the editing window.

Our vector system allows tandem integration of multiple **tRNA + target + gRNA** units (Fig. C), enhancing mutation efficiency. Plant endonucleases RNase P and RNase Z cleave the 5' and 3' ends of tRNA (Fig. A), releasing gRNAs that guide Cas9/gRNA to cleave the target DNA at the specified site.

In our lab, we modified the Tad-CBE vector by replacing the kanamycin resistance marker with hygromycin and incorporating the cotton endogenous U6 promoters (U6-7) for gRNA transcription. The protocol below focuses on constructing a vector with two target sites under the U6-7 promoter using rice Ubiquitin2 as the Cas9 promoter.

2. Target Search and Primer Design

CLA Sequence and Targets: Use the CRISPR-P website (<http://cbi.hzau.edu.cn/crispr/>) to search for targets (limited to *Gossypium raimondii* data). Alternatively, manually BLAST sequences for specificity. Prefer targets with unique restriction sites near the PAM (3 bp upstream) for downstream validation. G-to-A mutation on the negative strand corresponds to an C-to-T mutation on the positive strand. The target base is ideally located within positions 4 to 8 distal from the PAM site.

CTCGACGGACCTATACCACCTGTGGGAGCTTTGAGCAGTGCTCTCAGCAGGCTGCAAT

CAAACAGGCCTCTTAGAGAACTGAGAGAGGTTGCAAAGGTAAGTAGAAATGAATAGA
 ACCCAGAAAATAATATAATTTTACTAGAACTAAACTTACTAAGACTGGATGTGGT
 TTGATTGCAGGGAGTTA **CCAAGCAAATCGGTGGCCTATGCATGA**ACTGGCTG

CAAAAGTTGATGAGTATGCTCGTGGAATGATCAGTGGTTCTGGATCTACACTTTTCGAA
 GAACTTGGACTATATTATATTGGACCTGTTGATGGCCACAACATCGATGATTTAGTTTCT
 ATTCTCAAAGAGGTTAAGACTACTAAAACAACGGGTCCGGTCTTGATCCATGTTGTCAC
 TGAGAAAGGCCGAGGTTATCCATATGCAGAGAGAGCTGCTGACAAGTATCACGGTAAC
 ATACAGAATAAGCCTTTTTAGAAAGAGAGATCATTT **CATATGATTTAATTTAC** **TGG**TGCC

TCGATATCTGACCTTAGTTTAGGGAATAAAAGAATAAACTGTACCTGGATTTCTTTTCTC
 AGGAGTGGTGAAGTTCGATCCGGCAACTGGAAAGCAATTCAAAGGCAA
 TTCTGCTACCCAGTCTTACACTACATATTTTGCTGAGGCTTTGATTGCGGAAGCTGAGG
 CAGACAAAAATATTGTTGCCATCCATGCTGCAATGGGAGGTGGAACCGGATTAAACCT
 CTTCTCCGCCGGTCCCACAAAGATGTTTCGATGTGGGGATAGCTGAACAACATGCTG
 TCACCTTTGCTGCAGGCTTGGCCTGTGAAGGCTTGAAACCTTTTTGTGCAATCTACTCA
 TCATTCATGCAAAGGGCTTATGACCAG

Insert the target sequence and ligate between the MfeI digestion sites shown in the box of the vector diagram. The final ligation method is the one-step cloning method.

aagcatcagatgggca **AACAAAGCACCAGTGGTCTAGTGGTAGAATAGTACCCTGC**
CACGGTACAGACCCGGGTTTCGATTCCCGGCTGGTGCA *CATAGGCCACCGA*
TTTGCT **GTTTTAGAGCTAGAAATAGCAAGTTAAATAAGGCTAGTCCGTTA**
TCAACTTGAAAAAGTGGCACCGAGTCGGTGC **AACAAAGCACCAGTGGT**
CTAGTGGTAGAATAGTACCCTGCCACGGTACAGACCCGGGTTTCGATTCCCG
GCTGGTGCA *TTTCATATGATTTAATTTAC* *gtttagagctagaa*..... (载体 gRNA)

notes: Note: Lowercase = vector backbone; underscores = tRNA; bold = gRNA; italics = target site.

Primers::

CLA1 as: 5' AGCAAATCGGTGGGCCTATgtgcaccagccgggaat3' (靶标 1 互补序列)

CLA 2 s: 5' CATAGGCCACCGATTTGCTgttttagagctagaaata 3' (靶标 1 序列)

CLA 2 as: 5' GTAAATTAAATCATATGAAAgtgcaccagccgggaat3' (靶标 2 互补序列)

inf CLA as: 5' ttctagctctaaac GTAAATTAAATCATATGAAA3' (靶标 2 互补序列)

pRGEB32-7 s: 5' AAGCATCAGATGGGCAAACAAAGCACCAGTGGTCTAG3'

inf pRGEB32-7 s: 5' AAGCATCAGATGGGCAAACAAA 3'

The lowercase part is a fixed sequence. Primer design only requires modifying the target sequence indicated by the underline in the primer.


Among them, inf pRGEB32-7 s, pRGEB32-7 s, u6-7 s are universal primers and do not need to be changed. pRGEB32-7s contain a part of the vector sequence, facilitating In-fusion ligation. Because there are multiple trnas in the inserted sequence, inf pRGEB32-7 s, in order to enhance amplification specificity, have shorter sequence lengths than pRGEB32-7 s. Inf CLA as contains carrier sequence

connectors for convenient subsequent In-fusion connections.


3. First PCR Amplification

Overlap Extension PCR:

This method generates fusion fragments by designing primers with complementary overhangs. Two overlapping PCR products are combined and extended to form a single construct.

Overlapping extension schematic diagram

This is a schematic diagram of the test process

The gRNA+tRNA combination sequence can be synthesized by itself and ligated to a common vector (such as using the T vector). In the first PCR: Fragments 1 and 2 were amplified respectively from the bacterial liquid or plasmid containing the gRNA+tRNA vector. (The relevant sequence can be found in the attachment.)


Reaction Mixture (20 μ L):

	1	2
ddH ₂ O	16.1	16.1
Buffer	2	2
dNTP	0.3	0.3
S primer	pRGE32-7 s 0.2	CLA 2 s 0.2
AS primer	CLA1 as 0.2	CLA 2 as 0.2

Taq	0.2	0.2
Template DNA	1	1

PCR Conditions:

1. Initial denaturation: 95° C for 4 min
2. Denaturation: 95° C for 30 s
3. Annealing: 55° C for 30 s
4. Extension: 72° C for 20 s
 - Repeat steps 2 – 4 for 3 cycles
5. Denaturation: 95° C for 30 s
6. Annealing: 60° C for 30 s
7. Extension: 72° C for 20 s
 - Repeat steps 5 – 7 for 27 cycles
8. Final extension: 72° C for 5 min

Electrophoresis: Verify amplification.


4. Second PCR and Purification

Combine the two PCR products from Step 3 using overlap extension,

Reaction Mixture (100 μ L)		
ddH₂O		83.5μl
Buffer		10μl
dNTP		1.5μl
S primer	Inf pRGEB32-7 s	1μl
AS primer	Inf CLA as	1μl
Taq		1μl
Fragment 1		1μl
Fragment 2		1μl

PCR Conditions:

1. Initial denaturation: 95° C for 4 min
2. Denaturation: 95° C for 30 s
3. Annealing: 59° C for 30 s
4. Extension: 72° C for 20 s
 - Repeat steps 2 – 4 for 28 cycles
5. Final extension: 72° C for 5 min

Purification: Use a gel extraction kit to purify the product and measure concentration.

5. Vector Digestion

Digest the vector with BsaI.

Reaction Mixture (100 μ L):

	Volume (μL):
Tad-CBE	10 μ g
CutSmart Buffer	10 μ L
BsaI	1 μ L
ddH ₂ O	Up to 100 μ L

Digestion: 37°C for 5.5 h.

6. In-Fusion Cloning

Assemble the digested vector and insert using In-Fusion technology.

Components:

- Insert (100 ng)
- Linearized vector (100 ng)
- Exnase (0.5 μ L)
- CE Buffer (1 μ L)

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

7. Transformation

Transform competent cells with the ligation product.

Steps:

Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7s: TGTGCCACTCCAAAGACATCAG + Inf XX AS


Positive clones are sequenced (GC-rich regions).

Electrophoresis: Confirm positive clones.

Positive detection electrophoresis pattern:

Construct a successful carrier diagram:

8. Agrobacterium Electroporation

Electroporate the verified construct into *Agrobacterium*.

u6-7 s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

9. Validation

Genomic DNA Extraction: PCR amplify the target region, clone, and sequence.