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
——Using CLA as an Example


1. Basic Principles

After the target sequence is cleaved by the CRISPR system, the organism initiates its own repair mechanisms, leading to mutations and silencing of the target gene function (as shown below). The Cas9/gRNA complex can target any DNA site containing the motifs **5'-N₂₀-NGG-3'** or **5'-CCN-N₂₀-3'** (where NGG is the PAM required for Cas9 recognition). Cas9/gRNA precisely cleaves the target DNA double strand approximately 3 bp upstream of the PAM.


The dCas9, obtained through the targeted mutation of Cas9, retains only its targeting activity. When directed to the promoter sequence of a gene, its fused transcriptional activation domain can facilitate the targeted upregulation of that gene.

Our vector system allows tandem integration of multiple **tRNA + target + gRNA** units (Fig. C), enhancing mutation efficiency. Plant endonucleases RNase P and RNase Z cleave the 5' and 3' ends of tRNA (Fig. A), releasing gRNAs that guide Cas9/gRNA to cleave the target DNA at the specified site.

In our lab, we modified the pRGEB32 vector by replacing the kanamycin resistance marker with hygromycin and incorporating the cotton endogenous U6 promoters (U6-7, U6-9) for gRNA transcription. The protocol below focuses on constructing a vector with two target sites under the U6-7 promoter using rice Ubiquitin2 as the Cas9 promoter.

2. Target Search and Primer Design

CLA Sequence and Targets: Use the CRISPR-P website (<http://cbi.hzau.edu.cn/crispr/>) to search for targets (limited to *Gossypium raimondii* data). Alternatively, manually BLAST sequences for specificity. Prefer targets with unique restriction sites near the PAM (3 bp upstream) for downstream validation.

[-500]-[+5]

TTATATGTTTAGTAACATTTTTAGTTCAGTCTAACTAACTAATATAGTATAGTATCTAAATATTAATTTGATA
AGAATAAATTTAATAATTGCAATATTTAAATATAAATAATATCTAAATATTAATTAAGTAAAAATAAACTTAA
TAATTAATAATATCTAATACAAATTTTCATAATATTTTGACGATGAAAACATCTCAGGTAAGCCGAGTATACT
GTAACCATGTTTATTTTATATAATATCCAAACACCCCCACCCCCCAATTAACATTTCTAGCTGACCTTCGA
TACCAGCCCCGACCCCCAAGTTAAATAATCCATTTACATCAACAGCATAAAGTAACGTAAGCCCACAATG
ATCAGATTGAGAACGATAACAGACAACAACCTTCCCTCAAGAAATAGGAACAAAAGGAAAGAAAAAA
AAGGTACCTTCTTCTAACGCCGTAATCTATCCATTTGTTTTTCCATTGACGTGGTCCAATCAGAGAACG
ACCTGCCACT

Insert the target sequence and ligate between the two Bsa I digestion sites shown in the box of the vector diagram. The final ligation method is the one-step cloning method.

aagcatcagatgggcaAACAAAGCACCAGTGGTCTAGTGGTAGAATAGTACCCTGC
CACGGTACAGACCCGGGTTCGATTCCCGGCTGGTGCACATAGGCCACCGA
TTTGCTGTTTTAGAGCTAGAAATAGCAAGTTAAATAAGGCTAGTCCGTTA
TCAACTTGAAAAAGTGGCACCAGTCCGGTGCAACAAAGCACCAGTGGT
CTAGTGGTAGAATAGTACCCTGCCACGGTACAGACCCGGGTTCGATTCCCG
GCTGGTGCAGTGAAGTTCGATCCGGCAACgttttagagctagaa..... (载体 gRNA)

notes: Note: Lowercase = vector backbone; underscores = tRNA; bold = gRNA; italics = target site.

Primers::

CLA1 as: 5' AGCAAATCGGTGGGCCTATGtgccaccagccggaat3' (靶标 1 互补序列)
CLA 2 s: 5' CATAGGCCACCGATTGCTgttttagagctagaaata 3' (靶标 1 序列)
CLA 2 as: 5' GTTGCCGGATCGAACTTCACtgccaccagccggaat3' (靶标 2 互补序列)
inf CLA as: 5' ttctagctctaaaacGTTGCCGGATCGAACTTCAC3' (靶标 2 互补序列)
pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAAGCACCAGTGGTCTAG3'
inf pRGE32-7 s: 5' AAGCATCAGATGGGCAAACAAA 3'

(U6-9

pRGE32-9 s: 5' CAGCACATAACTGGCAAACAAAGCACCAGTGGTCTAG3'
inf pRGE32-9 s: 5' CAGCACATAACTGGCAAACAAA 3')


The lowercase part is a fixed sequence. Primer design only requires modifying the target sequence indicated by the underline in the primer.

Among them, inf pRGE32-7 s, pRGE32-7 s, pRGE32-9 s, inf pRGE32-9 s, u6-7 s, and u6-9 s are universal primers and do not need to be changed. pRGE32-7s and pRGE32-9s contain a part of the vector sequence, facilitating In-fusion ligation. Because there are multiple trnas in the inserted sequence, inf pRGE32-7 s and inf pRGE32-9 s, in order to enhance amplification specificity, have shorter sequence lengths than pRGE32-7 s and pRGE32-9 s. Inf CLA as contains carrier sequence connectors for convenient subsequent In-fusion connections.


3. First PCR Amplification

Overlap Extension PCR:

This method generates fusion fragments by designing primers with complementary overhangs. Two overlapping PCR products are combined and extended to form a single construct.

Overlapping extension schematic diagram

This is a schematic diagram of the test process


The gRNA+tRNA combination sequence can be synthesized by itself and ligated to a common vector (such as using the T vector). In the first PCR: Fragments 1 and 2 were amplified respectively from the bacterial liquid or plasmid containing the gRNA+tRNA vector. (The relevant sequence can be found in the attachment.)

Reaction Mixture (20 μ L):

	1	2
ddH₂O	16.1	16.1
Buffer	2	2
dNTP	0.3	0.3
S primer	pRGEB32-7 s 0.2	CLA 2 s 0.2
AS primer	CLA1 as 0.2	CLA 2 as 0.2
Taq	0.2	0.2
Template DNA	1	1

PCR Conditions:

1. Initial denaturation: 95°C for 4 min
2. Denaturation: 95°C for 30 s
3. Annealing: 55°C for 30 s
4. Extension: 72°C for 20 s
 - Repeat steps 2–4 for 3 cycles
5. Denaturation: 95°C for 30 s
6. Annealing: 60°C for 30 s
7. Extension: 72°C for 20 s
 - Repeat steps 5–7 for 27 cycles
8. Final extension: 72°C for 5 min

**Electrophoresis:** Verify amplification.


4. Second PCR and Purification

Combine the two PCR products from Step 3 using overlap extension,

Reaction Mixture (100 μ L)			
ddH ₂ O		83.5 μ l	
Buffer		10 μ l	
dNTP		1.5 μ l	
S primer	Inf pRGEB32-7 s	1 μ l	
AS primer	Inf CLA as		1 μ l
Taq		1 μ l	
Cas9 vector		2 μ l	

PCR Conditions:

1. Initial denaturation: 95°C for 4 min
2. Denaturation: 95°C for 30 s
3. Annealing: 59°C for 30 s
4. Extension: 72°C for 20 s
 - Repeat steps 2–4 for 28 cycles
5. Final extension: 72°C for 5 min

Purification: Use a gel extraction kit to purify the product and measure concentration.

5. Vector Digestion

Digest the vector with BsaI (Cas9) .

Reaction Mixture (100 μ L):

Volume (μ L):	
PCR Product	10 μ g
CutSmart Buffer	10 μ L

BSA	1 μ L
ddH ₂ O	Up to 100 μ L

Digestion: 37°C for 5.5 h.

6. In-Fusion Cloning

Assemble the digested vector and insert using In-Fusion technology.

Components:

- Insert (100 ng)
- Linearized vector (100 ng, Cas9)
- Exnase (0.5 μ L)
- CE Buffer (1 μ L)

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

7. Transformation

Transform competent cells with the ligation product.

Steps:

Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

u6-9s: GTCAAAAACACTATCCCACATTGCTAA + Inf XX AS


Positive clones are sequenced (GC-rich regions).

Electrophoresis: Confirm positive clones.

Positive detection electrophoresis pattern:

Construct a successful carrier diagram:

8. PCR amplification and purification of target bands:

Fragments containing the target sequence were ligated using PCR amplification with the enzymatic sites Hpa1 and Sbf1 at either end.

100μl 体系			
ddH ₂ O		83.5μl	
Buffer		10μl	
dNTP		1.5μl	
S primer	Inf TV s		1μl
AS primer	Inf sbf as		1μl
Taq		1μl	
Last vector		2μl	

Inf TV s : 5' CTGCGCACGTGTGTTaacctcgttttcacatgtg 3'

Inf sbf as : 5' ATTTGTGGACCTGCAgctaaaaacggactagcctt 3

PCR:

Initial denaturation: 95°C for 4 min

Denaturation: 95°C for 30 s

Annealing: 54°C for 30 s

Extension: 72°C for 40 s

Repeat steps 2–4 for 28 cycles


Final extension: 72°C for 5 min

Purification: Use a gel extraction kit to purify the product and measure concentration.

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10. In-Fusion Cloning

Assemble the digested vector and insert using In-Fusion technology.

Components:

- Insert (100 ng)
- Linearized vector (100 ng, dCas9-TV)
- Exnase (0.5 μ L)
- CE Buffer (1 μ L)

Incubation: 37°C for 30 min, then ice for 5 min. Store at -20°C.

11. Transformation

Transform competent cells with the ligation product.

Steps:

Mix with SOC medium and incubate at 37°C for 45 min.

Plate on kanamycin agar plates and incubate overnight.

Pick colonies for PCR validation using primers:

u6-7s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

u6-9s: GTCAAAAACCTATCCCACATTGCTAA + Inf XX AS

Positive clones are sequenced (GC-rich regions).

Electrophoresis: Confirm positive clones.

u6-7 s: TGTGCCACTCCAAAGACATCAG

Ubi-as: TGTTGGTCGCCGTTAGGA

12. Agrobacterium Electroporation

Electroporate the verified construct into *Agrobacterium*.

u6-7 s: TGTGCCACTCCAAAGACATCAG + Inf XX AS

Ubi-as: TGTTGGTCGCCGTTAGGA