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Mb2Cas12a protocol (英文版)

CRISPR/Mb2Cas12a Gene Editing Vector Construction (Kan) Method

(Taking CRISPR/Mb2Cas12a-tRNA as an example):

Digestion of Skeleton Vector Mb2

The skeleton vector Mb2 was digested with BsaI nuclease (NEB). The digestion effect was detected by electrophoresis, and the product was purified using a gel recovery kit followed by concentration measurement.

Table 1. Digestion System

100μl 体系	
Mb2	10μg
10X cut smart Buffer	10μl
BSA1	4μl
ddH ₂ O	up to 100μl

2. Single Target Vector Construction: tRNA-crRNA scaffold-gRNA-tRNA-TTTTTTT

- Method 1: Gene synthesis by Nanjing 金斯瑞 (Genscript).
- Method 2: Laboratory construction using Gh360 (kan) as the skeleton vector to amplify the tRNA-crRNA scaffold fragment.

2.1 First PCR Amplification: Amplify Target Sequence

Primers:

1041 F: AAGCATCAGATGGGCAAACAAAGCACCAAGTGGTCTAG

1041R:

CGAAAATAGCCATCGCCCCATGAATCTACCATAGTAGAAATTATGCACC

(Yellow: Target reverse complementary sequence; Green: crRNA scaffold; Red: tRNA)

The yellow part should be designed according to the user's target sequence.

PCR System:

ddH₂O	16.1	
Buffer	2	
dNTP	0.3	
S primer	1041F	0.2
AS primer	1041R	0.2
Taq	0.2	
模板	1	

PCR Conditions:

Pre-denaturation: 95°C, 5 min

Denaturation: 95°C, 30 s

Annealing: 55°C, 30 s

Extension: 72°C, 20 s

Cycles: 3 cycles

Denaturation: 95°C, 30 s

Annealing: 60°C, 30 s

Extension: 72°C, 20 s

Cycles: 27 cycles

Final extension: 72°C, 5 min

Electrophoresis detection required.

2.2 Second PCR Amplification: Add Recombinant Adaptor Sequence

Using the first PCR product as a template.

Primers:

1042 F: AAGCATCAGATGGGCAAACAAA

1042 R: CACTGGTGCTTTGTT**CGAAAATAGCCATCGCCCCATGA**

(Red: tRNA; Yellow: Target reverse complementary sequence)

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Cycles: 28 cycles

Final extension: 72°C, 5 min

Electrophoresis detection required.

In-fusion Ligation

Ligate the final target sequence into the BsaI site of the Mb2 vector.

Ligation System:

Component Volume (μl)

Target Fragment	100ng (2ul)
Linearized Expression Vector	100ng (1ul)
Exnase	0.5μl
CE Buffer	1μl
ddH ₂ O	Up to 5 ul

Incubate at 37°C for 30 min, place on ice for 5 min, and store at -20°C for later use.

Heat shock transformation, pick single colonies for sequencing, and preserve the correct plasmids in bacterial cultures.

3. Dual (Multiple) Targets: tRNA-crRNA scaffold-gRNA-tRNA-crRNA scaffold-gRNA-tRNA-TTTTTT

Based on the single-target vector, perform overlap extension PCR sequentially.

Sequences:

tRNA:

AACAAAGCACCAGTGGTCTAGTGGTAgaatagtaccctgccacggtacagacccgggTTCG
ATTCCCGGCTGGTGCA

crRNA scaffold:

TAATTTCTACTATGGTAGAT

Target: 23 bp

Attachments:

1. Mb2Cas12a skeleton vector (Mb2)
2. Gh360 template vector
3. Mb2Cas12a target design
4. Reference sequence for the final recombinant vector (Mb2Cas12a-Bsa1 after construction)